

BGZ-120 Auto Desiccant-supplying Machine



- (1).The machine is simple in structure, stable in function and widely applicable.
- (2).Apparatus available to stop the machine when no bottles detected. Safety control device available in case no desiccant is detected or desiccant jams.
- (3).All the processes are controlled by PLC program. The front operation panel of touching type is easy to operate, simple to maintain, convenient to clean and durable.

Specifications and Parameters

Max production capacity(Bottle/min): 50-60

Specification of Bottle(ml): 40-200

Voltage/Frequency(V/Hz): 220/50

Power(kW): 0.35

Air pressure (Mpa): 0.5-0.6

Air consumption(M³/min): 0.003

Measurement(mm): 1200×500×1500

Weight(kg): 200



Pharmaceutical Equipment Specialists

35 Bridgegate Drv, Victoria Dock, Hull, HU9 1SY, UK. tel: +44 1482 225628 fax : +44 8082089525

e mail: gorclark@gmail.com web: www.nevos.co.uk