

## DPH200 Aluminum-plastics blister Packing Machine



It apply to capsule,white pill,sugar-coat pill micelle,special-pill,marble chocolate,and chewing gum and so on. It is up to GMP standard.

Different fittings,it can be most satisfied with demand of users.

### Character

- Small volume,structure compact,high efficient
- Board-type pressure blow molding,roller-type sealing,beautiful production
- Have function of automation beating lot number tangent line
- Replacing mold and cutting are easy
- Non-side cutting,save materials,move smoothly,easy operation
- All of PLC automation control,frequency stepless speed regulation adjust range widely low-noise
- Close structure design,material use stainless steel being up to GMP standard

## Parameters:

max forming area	200×210mm	max width	210mm
max forming depth	12mm	inside diameter	φ70mm-φ76mm
cutting speed	Capsule 180 – 240ps/min	air consumption	50Liter/min
	tablet 240 – 300 ps/min	water consumption	3.2Liter/min
cutting distance	105mm	Sum-efficient	8.0KW
forming distance	200mm	power source	380V 50HZ
max outside diameter	PVC φ400mm	Weoght	1500KG
	PTP φ280mm	size kf contour (length width height)	2400×900×1650mm



*Pharmaceutical Equipment Specialists*

35 Bridgegate Drv, Victoria Dock, Hull, HU9 1SY, UK. tel: +44 1482 225628 fax : +44 8082089525

e mail: [gorclark@gmail.com](mailto:gorclark@gmail.com) web: [www.nevos.co.uk](http://www.nevos.co.uk)