

# PZY2000I Tablet and capsule weight checker

The equipment is mainly used for pharmaceutical industry and inspector of pharmaceutical product line.



## Main characteristic:

1. The equipment has favorable man machine conversation interface to set parameters.  
Parameter setting: demarcate of weighing table; sampling time and interval time; sending medicine time; standard weight setting of single tablet or capsule; print time; code set and change.
2. Parameter can be changed according to produce demand at any moment.
3. The equipment can sampling, weigh, contrast and analyze data, draw dynamic curve and form, print and query automatically.
4. Its construction is compact, testing is steady and precise, and minor error.
5. The equipment has absorbing-dust device which can absorb dust of tablet and capsule before weighing to ensure checking accurately.

## Main technical data:

Max weight of capsule or tablet: 350 g

Error: 0.5%

Power supply: AC220V50Hz

Power: 0.2KW

Machine size (W×D×H) : 1000×500×1100(mm)



*Pharmaceutical Equipment Specialists*

35 Bridgegate Drv, Victoria Dock, Hull, HU9 1SY, UK. tel: +44 1482 225628 fax : +44 8082089525

e mail: [gorclark@gmail.com](mailto:gorclark@gmail.com) web: [www.nevos.co.uk](http://www.nevos.co.uk)

